INSTRUCTIONS FOR USE / INSTRUCCIONES DE USO

Active HEX

EN (ENGLISH).....	2
ES (ESPAÑOL)	15
FR (FRANÇAIS)	29
DE (DEUTSCH)	42
IT (ITALIANO).....	56
PT (PORTUGUÊS).....	70
GR (Ελληνικά).....	83
NL (NEDERLANDS)	97
PL (POLSKI)	111
CZ (ČEŠTINA)	124
HR (HRVATSKI)	136
HU (MAGYAR)	148
BG (БЪЛГАРСКИ)	161
LT (LIETUVIŲ K.).....	175
LV (LATVIEŠU).....	188
SK (SLOVENČINA)	202
RO (ROMÂNĂ).....	216
SR (SRPSKI)	231
DA (DANSK)	244
FI (SUOMI).....	258
NO (NORSK).....	271
SV (SVENSKA)	283
ET (EESTI).....	296
(العربية) AR	308
SL (SLOVENŠČINA).....	321

EN (ENGLISH)

Active HEX implants

Rx Only

Caution: Federal law in the USA restricts the sale of this device to licensed healthcare professionals.

Please read these instructions carefully before using the implants and components supplied by Terrats Medical S.L. The following instructions are intended for qualified professionals who are specialists in the field of dental implantology and have received the appropriate training.

The user must ensure that the selected product is suitable for the intended purposes and procedures.

Instructions for use are also available on the following website ifu.dessdental.com/activeheximplant

Description of the product

The DESS® Active HEX dental implant is part of an endo-osseous subgingival connection implant system consisting of dental implants, prosthetic components and surgical instruments. DESS® implants are made of titanium (grade 4) according to the ISO 5832-2 standard, with a surface topography obtained by a surface treatment consisting of sandblasting and double acid etching. DESS® Active HEX implants are available in a range of diameters from 3.0 – 5.5.

Intended use and indications

The product is intended to be inserted surgically into the bone of the upper or lower jaw to support other prosthetic products, such as an artificial tooth, to restore the mastication function of the patient.

DESS® implants are intended for use in the upper or lower jaw in patients who are partially or completely edentulous. They are designed to support single and multi-unit restorations as well as to retain overdentures. Dental implants can be used for immediate or delayed loading techniques. Immediate loading is only recommended when adequate primary stability is achieved. The smallest diameter implants, from 3.0 to 3.5 mm, are intended for use in reduced interdental spaces, where there is not enough alveolar bone for an implant with a larger diameter. The use of small diameter implants in posterior rehabilitation is not recommended. Larger diameter implants are intended for both upper and lower jaw rehabilitations for functional and aesthetic rehabilitations in partially or totally edentulous patients.

Introducing the DESS® dental implant

DESS® dental implants undergo a meticulous manufacturing, control and cleaning process before being packaged in a sterile area and subsequently sterilised by irradiation. The implants are supplied in a sealed vial, inside which the implant is secured using titanium mounts to avoid contact with any material other than titanium and to avoid possible contamination of its surfaces. The vial is placed in a heat-sealed blister to keep it sterile until its expiry date. Both the outer package (box) and the blister have a label with the following information: the batch number, size and model of the implant and its expiry date.

To ensure the mandatory traceability to which the dental implant is subject, retain the label or transcribe its data to the patient's file.

DESS® dental implants are supplied sterile and are intended for single use only. The implants must not therefore be re-sterilised or reused under any circumstances.

Terrats Medical S.L. will accept no responsibility for re-sterilised implants, regardless of who performed the re-sterilisation or what method was used.

Contra-indications

DESS® dental implants should not be placed in patients who are found to be medically unsuitable for the intended treatment.

Contra-indications include, but are not limited to: titanium allergies, metabolic or systemic disorders associated with wound and/or bone healing, previously irritated bones, use of pharmaceuticals that inhibit or alter natural bone remodelling, uncontrolled diabetes, alcoholism or drug addiction, blood clotting disorders, anticoagulant therapy, diseases with periodic use of high doses of steroids, metabolic bone disease, chemotherapy or radiotherapy, chronic periodontal inflammation, insufficient soft tissue coverage, any disorder that inhibits the patient's ability to maintain adequate daily oral hygiene, uncontrolled parafunctional habits, moderate or severe smoking, uncontrolled endocrine

diseases, insufficient bone height and/or width, insufficient interarticular space, pregnancy, psychosis and vascular conditions.

DESS dental implants are not recommended for children under 18 years of age. Treatment of youngsters (of 18 to 21 years of age) is not recommended until skeletal growth is complete and epiphyseal closure has occurred. The use of unloading splints is recommended in patients with parafunctional habits.

Intended users and patient groups

DESS® implantology products should only be used by dental specialists with experience in maxillary implantology and other specialties, such as dental diagnostics, planning, dental surgery or prosthetic techniques. Their use is limited to dental laboratories and clinics.

Indicated for edentulous patients (totally or partially) who require oral rehabilitation by means of implant-supported prostheses. Rehabilitation can be unitary, multiple or overdentures, both in the upper and lower jaw. Its use is indicated in patients with complete dentofacial development. There are no differences in end users by age group, sex, ethnicity, family predisposition or genetic aspects. It is not indicated for use in patients without dental problems.

Clinical benefits and side effects

As a clinical benefit, patients can expect total or partial rehabilitation of the dentition, allowing them to regain proper masticatory function.

It is the surgeon's responsibility to provide the patient with any and all information about side effects, precautions and possible complications that may arise after implant surgery, as well as to complete the consent form.

Warnings and precautions

DESS® dental implants should be stored in their original protective packaging, in a cool and dry place (at room temperature) protected from direct exposure to sunlight.

DESS® dental implants must not be reused under any circumstances. Reuse of products labelled single-use may also result in the device not functioning as intended due to a change in its geometry. No warranty claim arising from the reuse of a single-use device will be accepted.

The use of electrosurgery is not recommended due to the conductivity of dental implants.

DESS® implants have not been evaluated for safety, compatibility, heating or migration in MRI environments. Prior to surgery, it is the physician's responsibility to check the condition of the implant packaging and whether it is the correct product to meet the patient's needs. Terrats Medical S.L. recommends always having additional implants available. When using our products during intraoral procedures, mandatory precautions must be taken to avoid accidental aspiration/ingestion by the patient.

As part of the dental implant rehabilitation and prior to surgery, the physician is solely responsible for conducting complete preoperative evaluation and planning of the entire treatment with the patient. Precise planning is necessary to avoid any accidental damage during bed preparation and implant placement to vital anatomical structures, such as the mental nerve and inferior alveolar nerve, as this may result in anaesthesia, paraesthesia or dysaesthesia.

Small diameter implants and angled abutments are not recommended for the posterior region

Serious Incident Warning

For patients / users / third parties in the European Union with an identical regulatory regime (Regulation 2017/745/EU): If any serious incident should occur due to the use of the product, notify the manufacturer and your national authority.

The manufacturer's contact information for reporting this type of incident is:

Terrats Medical SL
Calle Mogoda 75-99
08210 Barberá del Vallés, Barcelona (Spain)
Tel: (+34) 935646006
medical.incidents@dessdental.com

Surgical protocol

It is important to examine and evaluate each patient carefully to determine each particular situation, especially any bone or soft tissue deficits that may influence the final outcome. To create optimal healing conditions for the implant, the hard and soft tissues must be treated in a non-aggressive manner. The implant placement site must be carefully

prepared, ensuring the elimination of any contaminants or sources of infection. The following measures are recommended to minimise excessive temperatures.

Recommended drills:

Before starting the surgery, make a base point on the bone using the Ø1.5" diameter drill.

Platform	Soft bone (IV)	Medium bone (II-III)	Dense bone (I)
3.0	1.5	2.0	2.0 2.4/2.8
3.5 NP	2.0 2.4/2.8*	2.0 2.4/2.8 (2.8/3.2)*	2.0 2.4/2.8 2.8/3.2
4.3 RP	2.0 2.4/2.8 (2.8/3.2)*	2.0 2.4/2.8 3.2/3.6	2.0 2.4/2.8 3.2/3.6 (3.8/4.2)*
5.0 RP	2.0 2.4/2.8 3.2/3.6	2.0 2.4/2.8 3.2/3.6 3.8/4.2	2.0 2.4/2.8 3.2/3.6 3.8/4.2 (4.2/4.6)*
5.5WP	2.0 2.4/2.8 3.2/3.6 (3.8/4.2) *	2.0 2.4/2.8 3.2/3.6 3.8/4.2 4.2/4.6 (4.2/5.0) *	2.0 2.4/2.8 3.2/3.6 3.8/4.2 4.2/4.6 4.2/5.0 Screw Tap

*Dense cortex status

Drilling should be performed with alternating movements at a maximum speed of 2000 rpm without exceeding 70 Ncm, and constant irrigation is necessary to avoid bone overheating.

Bone conformers (screw taps) are intended for use in the upper or lower jaw to prepare the osteotomy in situations involving dense bone to avoid excessive torque during implant insertion (max 70 Ncm for NP, RP and WP and 45 Ncm for 3.0 implants).

Note. The actual implant length is 0.5 mm shorter than the drill mark.

Caution. The drills for DESS® Active HEX implants do not have 8.5 and 11.5 depth marks.

Recommended speeds for drills:

Diameter	Ø 1.5 m m	Ø 2.0 m m	Ø 2.4/2.8 mm	Ø 2.8/3.2 mm	Ø 3.2/3.6 m m	Ø 3.8/4.2 mm	Ø 4.2/4.6 mm	Ø 4.2/5.0 mm	Screw tap (3.0/3.5/4.3/5.0/5.5)
Speed	2000 Rpm	2000 Rpm	2000 Rpm	2000 Rpm	2000 Rpm	2000 Rpm	2000 Rpm	2000 Rpm	25 Rpm

Use sharp, sterilised instruments

Maximum number of recommended uses: 20

Maintain a low temperature

For drilling, an intermittent drilling technique is recommended, with high levels of refrigeration, using pre-cooled physiological saline solution or at room temperature.

Implant placement

1. Bone bed preparation. If the surgical protocol used is flapless, add the soft tissue height to the depth of the reaming.

2. Use the depth gauge to check the milling depth.
3. Remove the vial from the inside of the container.
4. Insert the implant insertion key previously installed in the handpiece directly into the dental implant connection.
5. Place the implant in the osteotomy previously performed and proceed with the insertion. Implant placement should be performed at very low speed (25-30 rpm).
It is recommended that the alveolus be full of blood clot at the moment of insertion, in order to impregnate the implant with blood cells.
6. For the placement of the implants never exceed the insertion torque of 45 Ncm in 3.0 implants or 70 Ncm for 3.5, 4.3, 5.0 and 5.5 implants.
Caution: The application of a higher insertion torque may damage the internal connection of the implant, and cause fracture or necrosis in the bone bed.
Implants can be tilted up to 45°, but any angulation greater than 30° should be splinted.
Great care and attention must be taken to prevent the patient from swallowing or aspirating the components and/or tools used.
After implant placement, the surgeon must evaluate bone stability and quality to determine when to load the implant.
7. Depending on the selected surgical protocol, place a cover screw, healing abutment, temporary restoration or Multiunit abutment with corresponding cover cap.

Storage, handling and transport.

The device should be stored and transported dry in its original packaging at room temperature and should not be exposed to direct sunlight. Incorrect storage or transport may affect the characteristics of the product.

Disposal

Disposal of the devices should follow local regulations and environmental requirements, taking into account different levels of contamination. Consult Terrats Medical for the implementation of product packaging material recycling.

Exclusion of warranties and liability

DESS® dental implants are part of a comprehensive concept and should only be used with original components and tools supplied by DESS®. The use of products manufactured by third parties which are not supplied for the DESS® implant will automatically void any express or implied warranty or liability.

Magnetic resonance imaging (MRI) safety information

Warning: The RF safety of the device has not been tested. A patient with this device can be safely scanned on an MRI system under the following conditions:

- Static magnetic field strength (B_0) less than or equal to 3.0T.
- Special gradient magnetic fields of maximum 3000 gauss/cm (30 T/m).
- RF excitation: Circular Polarization (CP)
- For the body coil, a reference point should be marked at least 30 cm from the implant or ensure that the implant is outside the coil. Extremity T/R coils are allowed. Head T/R coil excluded.
- Normal mode of operation in the allowed imaging zone.
- Maximum whole body specific absorption rate of 2 W/kg (Normal mode of operation).
- Maximum head specific absorption rate not evaluated.
- Scan duration without specific limitations due to implant heating.

Note: removable restorations must be removed prior to scanning.

Summary of Safety and Clinical Performance (SSCP)

The Summary of Safety and Clinical Performance (SSCP) for products belonging to the families of devices is available in the European Database on Medical Devices (EUDAMED). You can consult the document linked to the basic UDI-DI 8435457217ST01P9Y on the EUDAMED public website: <https://ec.europa.eu/tools/eudamed/#/screen/home>

Active HEX surgical drills

OBJECT.

The objective of the surgical drills is to allow an implant alveolus to be made in the patient's bone support. The realization of the implant alveolus will be done by the sequential step of drills of different diameters.

DESCRIPTION.

The surgical drills are made of stainless steel suitable for this use. The drills are supplied non-sterile. They can be used several times; it is recommended not to exceed 20 drilling cycles and they should be thoroughly cleaned and sterilized before each use as indicated in point 4.

The burs have a conical shape at the cutting end and a contra-angle handpiece connection at the other end.

The burs have markings to control the drilling depth. Drilling should be carried out with alternate movements, not exceeding 70 Ncm, and constant irrigation is necessary so that the bone does not heat up and can cause bone necrosis. The drilling sequence described in point 3 of the drilling protocol must be followed.

The drill must not be stopped while drilling is in progress, as it could become blocked and break during extraction.

The drills do not require special storage conditions but should be kept in a clean and dry place.

Burs that are no longer in use must be cleaned to avoid pathogenic residues, destroyed and treated as non-hazardous ferrous waste.

Adverse effects

If the drilling sequence is not respected, if adequate irrigation is not carried out, if cleaning and sterility are not carried out correctly, implant procedure may fail due to lack of osseointegration of the implant, and bone necrosis or serious oral infections may occur.

Precautions

Do not leave the burs loose in the metal box to avoid knocks.

Do not let residues (blood, tissue, etc.) dry on the instruments. Clean immediately after surgery. See section 4 in this document for more details.

Do not store wet instruments in the surgical box.

SURGICAL PROTOCOL

Depending on the density of the bone (type 1 = very hard bone, type 4 = very soft bone) different drilling protocols should be applied for the Active HEX implant. This will make the adjustment of the implant in the alveolar preparation to the individual bone types more flexible depending on their quality and anatomical situation.

Drilling sequence

The burs are made of stainless steel with a superficial DLC (Diamond like Carbon) finish. The burs are used with external irrigation and are available in 3 lengths:

- 7-10 mm (28 mm)
- 7-15 mm (33 mm)
- 10-18 mm (36 mm)

Use an insertion and extraction movement and drill the bone for 1-2 seconds. Pull the drill upwards without stooping the surgical motor. This allows irrigation to remove possible bone chips.

Proceed until the desired drilling depth is reached.

Screw taps are available for dense bone situations to avoid excessive torque during implant insertion (max 70 N-cm for NP, RP and WP and 45 N-cm for 3.0 implants). Do not exceed the maximum speed of 2000 rpm.

Considerations for soft bone

The self-grinding capability of the Active HEX implant allows it to be inserted in places that have been prepared with little depth. This capability gives the implant a lot of usability in situations where we have nearby vital anatomical structures or in soft bone when maximum bone consideration is desirable. Drill to 2-4 mm less than the total length, insert the implant at the drilled depth and continue with insertion. The implant will drill its own way to the final depth.

Considerations for dense bone

Self-milling with the implant itself should not be done with dense bone. The tapping should be used for dense bone when the standard drilling protrusion is not sufficient for full implant insertion without exceeding the maximum recommended insertion torque (max. 70 Ncm for NP, RP and WP) and 45 Ncm for 3.0 implants. Do not exceed the maximum speed of 2000 rpm.

DESS Active HEX 3.0 mm

Recommended steps

Dense cortex only, not drill to the full drilling depth

Diameter implant	Type of bone	Diameter drill			
		Ø 1.5	Ø 2.0	Ø 2.4 / 2.8	Screw tap
3.0	Very hard bone I	-			
	Hard bone II	-		-	-
	Soft bone III	-		-	-
	Very soft bone IV		-	-	-
	RPM Max	2000	2000	2000	25

DESS Active HEX 3.5 mm NP

Recommended steps

Dense cortex only, not drill to the full drilling depth

Diameter implant	Type of bone	Diameter drill			
		Ø 2.0	Ø 2.4 / 2.8	Ø 2.8 / 3.2	Screw tap
3.5	Very hard bone I				
	Hard bone II				-
	Soft bone III				-
	Very soft bone IV			-	-
	RPM Max	2000	2000	2000	25

DESS Active HEX 4.3 mm RP

Recommended steps

Dense cortex only, not drill to the full drilling depth

Diameter implant	Type of bone	Diameter drill					
		Ø 2.0	Ø 2.4 / 2.8	Ø 2.8 / 3.2	Ø 3.2 / 3.6	Ø 3.8 / 4.2	Screw tap
4.3	Very hard bone I			-			
	Hard bone II			-		-	-
	Soft bone III			-		-	-
	Very soft bone IV				-	-	-
	RPM Max	2000	2000	2000	2000	2000	25

DESS Active HEX 5.0 mm RP

Recommended steps

Dense cortex only, not drill to the full drilling depth

Diameter implant	Type of bone	Diameter drill					
		Ø 2.0	Ø 2.4 / 2.8	Ø 3.2 / 3.6	Ø 3.8 / 4.2	Ø 4.2 / 4.6	Screw tap
5.0	Very hard bone I						
	Hard bone II					-	-
	Soft bone III					-	-
	Very soft bone IV				-	-	-
	RPM Max	2000	2000	2000	2000	2000	25


DESS Active HEX 5.5 mm WP

Recommended steps

Dense cortex only, not drill to the full drilling depth

Diameter implant	Type of bone	Diameter drill						
		Ø 2.0	Ø 2.4 / 2.8	Ø 3.2 / 3.6	Ø 3.8 / 4.2	Ø 4.2 / 4.6	Ø 4.2 / 5.0	Screw tap
5.5	Very hard bone I					-		
	Hard bone II							-
	Soft bone III							-
	Very soft bone IV					-	-	-
	RPM Max	2000	2000	2000	2000	2000	2000	25

CLEANING, DISINFECTION AND STERILIZATION

Cleaning.

The instruments should be cleaned manually after use. To do this, brush and rinse with water the excess dirt and particles on the entire surface of the instruments for 25 to 35 seconds.

Disinfection.

Immerse the instruments in a disinfectant bath suitable for dental material. Strictly follow the instructions for use and prescriptions of the selected disinfectant. We recommend the disinfectant specially indicated for dental sanitary equipment Instrument (universal disinfectant for instruments).

Rinse thoroughly with water after application of the disinfectant.

Warning.

Do not use solutions containing ammonia, hydrogen peroxide or acidic substances, as these can damage the surface coating of the burs.

Sterilization of the units.

Sterilize the metal product in a steam autoclave. It is recommended that a sterilization cycle of 134°C and a minimum of 6 minutes be used. Wait until the drying cycle is completely finished. It is recommended that sterilization controls be used, recording the date and expiration date, in addition to carrying out periodic controls of the sterilization process using biological indicators.

The presence of corrosion after sterilization is the main factor in rejecting the use of instruments, independently of whether or not they have cutting capacity. Inspect the instruments once they have been sterilized in order to verify if there is depreciation after the sterilization cycles.

The burs should not be sterilized in their original packaging; use the specific bags for sterilization.

GLOSSARY OF SYMBOLS

Batch

Reference

Unique device identifier

Model number

Medical Device

Sterilised by radiation

Single sterile barrier system with the protective packaging inside

Expiry date

Do not reuse

Do not use if the packaging is damaged

Keep out of sunlight

Non sterile

Do not re-sterilize

CE marking with Notified Body number

Caution

Date of manufacture

Rx ONLY

Prescription only

MR Conditional

Please refer to the Instructions for Use
Download link
<https://www.dessdental.com/en/downloads>

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UKRP

Cyberpan Support Services Ltd
 The UnionSuite,
 The Union Building,
 51-59 Rose Lane, Norwich, Norfolk
 England NR1 1BY

IMPLANT CARD

DESS® implants include an implant/patient card to be filled in by the practitioner and given to the end patient. This card allows the traceability of the implants used for treatment to be maintained.

Steps to follow:

1. Fill in the front side of the card with the patient's details, the date of surgery and your own professional details.
2. Unhook the main label from the blister pack and attach it to the back of the card.

	Title	Description
	Patient identification	Indicates the patient's personal information
	Patient information website	Indicates a website where a patient can obtain additional information on the medical product
	Health care center or doctor	Indicates the Health Care Centre/doctor's address where medical information about the patient may be found
	Date	Indicates the date the information was entered or the Medical procedure took place

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